




ORDERING GUIDE

REUSABLE

CODE	DESCRIPTION	N°/PACK
EU5086	ECMOLIFE SYSTEM _____	1
	includes: Console _____	1 (A)
	Motor driver _____	2 (B)
	Flow meter 3/8" x 3/32" _____	2
	Bracket _____	1 (C)
	Venous probe _____	1 (D)
	Bubble sensor 3/8" x 3/32" _____	1
	ECMOLife power cord EU _____	2
EU5093	LEONARDO Trolley _____	1 (E)
	includes: ECMOLIFE Cable Holder _____	1 (F)
	REOMED 1000 isolating trans-formers 1000VA _____	1
	REOMED Plug C13-C14 _____	2
EU5094	LEONARDO Slim _____	1
	includes: ECMOLife Cable Holder _____	1 (F)
	REOMED 1000 isolating trans-formers 1000VA _____	1
	REOMED Plug C13-14 _____	2



REUSABLE		
CODE	DESCRIPTION	N°/PACK
R229P01-010	MAX BLEND 0-15 LPM	1 *
R100P26	Bracket, Pole mount, gas blender	1 *
R129P91	Air Tube	1 *
R129P94	Oxygen Tube	1 *
R219P31	3 in 1 single pack (connector gas outlet)	1 *
EU3941	ECMOLIFE HEATER COOLER	1
	includes: HC Unit _____	1
	UK AND EU Power Supply Cables _____	1
	Hoses 1,5m _____	2
	Hanses connection tube _____	1
	Water Discharge Hose (Female) _____	1
EU3943	Hoses 3m	2
EU3898	Ambulance Plate holder	1
EU3996	FOIL Transport System	1
EU1900	ECMOLIFE FLOWMETER 3/8" x 3/32" (only for ECMOLIFE 1.0)	1
EU6726	ECMOLIFE FLOWMETER/BUBBLE SENSOR 3/8" x 3/32" (only for ECMOLIFE 2.0)	1
EU6727	ECMOLIFE BUBBLE SENSOR 3/8" x 3/32" (only for ECMOLIFE 2.0)	1

REUSABLE		
CODE	DESCRIPTION	N°/PACK
EU6728	ECMOLIFE FLOWMETER 1/4" x 3/32" (only for ECMOLIFE 2.0)	1
EU6729	ECMOLIFE FLOWMETER/BUBBLE SENSOR 1/4" x 3/32" (only for ECMOLIFE 2.0)	1
EU1901	ECMOLIFE FLOWMETER 1/4" x 3/32" (only for ECMOLIFE 1.0)	1
EU3953	AMBULANCE 24V POWER	1
EU6734	EDWARDS ECMOLIFE CABLE	1
EU10775	ECOMLIFE Power Cord GB	1 *
EU10996	ECMOLIFE Power Cord BR	1 *
2M85274	Holder with BAR 25x10 mm	1  *
EU3899	ECMOLIFE Auxiliary batteries Kit with transport bag	1 
EU2511	ECMOLIFE Holder Motor Driver	1 

* Product not manufactured by EUROSETS S.r.l.

• Products in bold are manufactured and certified by Eurosets S.r.l.

• ECMOLife System, ECMOLife Heater Cooler and Single-use codes are CE₀₁₂₃ Certified

SINGLE USE

CODE	DESCRIPTION	N°/PACK
30AG70232	ADULT ECMO tubing set Includes: adult oxy, adult/paed pump, circuit pressure transducer Metallux 3 to 1	1
30AG70233	PAEDIATRIC ECMO tubing set Includes: paediatric oxy, adult/paed pump, circuit pressure transducer Metallux 3 to 1	1
30AG70234	NEW BORN ECMO tubing set Includes: new born oxy, new born pump, circuit pressure transducer Metallux 3 to 1	1
30AG70232/SM	ADULT ECMO tubing set Includes: adult oxy, adult/paed pump, circuit pressure transducer Smith Medical	1
30AG70231/SM	VAS tubing set Includes: adult/paed pump, circuit pressure transducer Smith Medical	1
14AG70231/SM	ECMOLIFE VAS tubing set Includes: pump, circuit pressure transducer Smith Medical	1
14AG70232	Adult ECMO tubing set Includes: adult oxy, adult pump, circuit pressure transducer Metallux 3 to 1	1
14AG70232/SM	ECMOLIFE adult ECMO tubing set Includes: adult oxy, adult pump, circuit pressure transducer Smith Medical	1
14AG70233	Pediatric ECMO tubing set Includes: paediatric oxy, adult/paed pump, circuit pressure transducer Metallux 3 to 1	1
14AG70234	New born ECMO tubing set Includes: new born oxy, new born pump circuit pressure transducer Metallux 3 to 1	1
AG70232/PP	ECMOLife 6 hrs Tubing Set Includes: PP Oxy, pump, transducer Smith Medical	1
EU5083	ECMO ADULT console Sterile (Stand alone ECMO Oxy)	2
EU5084	ECMO PAEDIATRIC console Sterile (Stand alone ECMO Oxy)	2
EU5165	ECMO NEW BORN console Sterile (Stand alone ECMO Oxy)	2
AG5078	ADULT/PAEDIATRIC centrifugal pump Sterile (Stand alone)	1
AG5087	NEW BORN centrifugal pump Sterile (Stand alone)	1

**30
DAYS
VALIDATED**

**14
DAYS
VALIDATED**

**30
DAYS
VALIDATED**