


PHOSPHORYLCHOLINE COATING

Phosphorylcoline or PC is the predominant lipid headgroup found within the outer layer of cell membranes. PC has a natural affinity for water and binds water tightly around itself. As a result, the outer layer of the cell membrane does not promote clot formation (thrombosis).



ORDERING GUIDE

CODE	DESCRIPTION	N°/PACK
AG5032	TRILLY FULL COATED STERILE	2
EU2331	Holder for TRILLY	1



Prodotto da Manufactured by Fabriquè par Hersteller Fabricado por

Eurosets s.r.l. Strada Statale 12, n°143 41036 Medolla (MO) Italy Tel: +39 0535 660311 Fax: +39 0535 51248 E-mail: info@eurosets.com eurosets.com







