


Phosphorylcholine or PC is the predominant lipid headgroup found within the outer layer of cell membranes. PC has a natural affinity for water and binds water tightly around itself. As a result, the outer layer of the cell membrane does not promote clot formation (thrombosis).



ORDERING GUIDE

CODE	DESCRIPTION	NO./PACK
AG5232	TRILLY PEDIATRIC AF (Cardiotomy reservoir pre-connected with oxygenator module, PC coated)	2
EU2331	TRILLY PEDIATRIC AF Holder	1

REFERENCES

1. Eurosets performance data from lab tests.
2. Eurosets performance data from lab testing according to ISO-7199 conditions.

3. Richmond ME, Charette K, Chen JM, et al. The effect of cardiopulmonary bypass prime volume on the need for blood transfusion after pediatric cardiac surgery. *J Thorac Cardiovasc Surg* 2013;145:1058-64. 10.1016/j.jtcvs.2012.07.016 [PubMed] [CrossRef] [Google Scholar]
4. Lou S, Ji B, Liu J, Yu K, Long C. Generation, detection and prevention of GME during CPB. *Int J Artif Organs*. 2011 Nov; 34(11):1039-51
5. G.J. Meyers RRT, CCP. Does turning on the vent increase microemboli transmission post arterial filter? EBCC 2009 Wien

EUROSETS EUROPE

ITALY

Eurosets S.r.l.
Strada Statale 12, 143
41036 Medolla (MO)
Italy
Ph: +39 0535 660311
Fax: +39 0535 51248
info@eurosets.com

BENELUX

Eurosets Benelux
Winstar Park,
Rue Provinciale, 62
1301 Bierges (Wavre)
Belgium
Ph: +32 10237993
info.nl@eurosets.com
info.be@eurosets.com

FRANCE

Eurosets France Sarl
Immeuble LE FONTENAY
96, Boulevard Vivier
Merle
69423 LYON Cedex 03
Ph: +33 962191542
Fax: +33 970066573
info.fr@eurosets.com

GERMANY

Eurosets GmbH
Industriestr. 29
82194 Gröbenzell
Ph: +49 81426671227
Fax: +49 81426671229
info.de@eurosets.com

EUROSETS ASIA

CHINA

Eurosets (Suzhou) Medical Device Ltd.
Suite 1308-2, The Summit (South Block)
No. 119, Suzhou Avenue West,
Suzhou Industrial Park - 215000, Greater China.
苏州工业园区苏州大道西119
号苏悦广场南楼1308-2, 215000,
更大的中国
info.cn@eurosets.com

   eurosets.com

STERILE/EO



CE 0123

Latex free



Prodotto da
Manufactured by
Fabriqué par
Hersteller
Fabricado por

Eurosets s.r.l.
Strada Statale 12, n°143
41036 Medolla (MO), Italy
Tel: +39 0535 660311
Fax: +39 0535 51248
E-mail: info@eurosets.com
eurosets.com

Distribuito da
Distributed by
Distribué par
Vertrieb
Distribuido por